



<16.7 km
Range per Charge



136 kg
Weight Capacity



94 cm
Turning Radius



GO GO ELITE TRAVELLER[®] 2

Navigate Tight Corners with Ease!

The Go Go Elite Traveller 2 boasts the tightest turning radius of all 4-wheel scooters, measuring just 94 cm. With Pride's patented EZ Turn Technology, maneuvering through tight spaces has never been easier, enabling you to execute U-turns on sidewalks, in grocery store aisles, and even in elevators. Additionally, the scooter is equipped with a turning speed sensor to prevent turns at excessive speeds, ensuring safe and controlled maneuvering.



GOGO ELITE TRAVELLER[®] 2



Pride's EZ Turn Technology



Feather Touch Disassembly



Breathable Mesh Seat

RANGE PER CHARGE (MAX) ^{1,3}	Up to 16.7 km
MAXIMUM WEIGHT CAPACITY	136 kg
TURNING RADIUS ²	94 cm
MAXIMUM SPEED ¹	6.8 km/h
LENGTH & WIDTH ^{2,5}	L: 103.6 cm W: 49.4 cm
SEAT TYPE WIDTH & DEPTH (STANDARD)	Compact, foldable ⁸ W: 43 cm (17") D: 43 cm (17")
SEAT-TO-DECK HEIGHT ²	41.6 cm - 45.3 cm
SEAT-TO-GROUND HEIGHT ²	46.3 cm - 51.4 cm
GROUND CLEARANCE ^{2,7}	3.4 cm at transaxle

OVERALL WEIGHT (NOT INC BATTERY) ^{2,6}	37.6 kg
WEIGHT OF HEAVIEST PIECE ²	15.8 kg (front section)
SUSPENSION	N/A
3W/4W OPTION	4W
FRONT TYRES	2" x 7" solid
REAR TYRES	2.5" x 8" solid
BATTERY REQUIREMENTS ⁹	(2) 18AH
BATTERY CHARGER	Off-board, 2A
BATTERY WEIGHT ²	6 kg
COLOUR OPTIONS	■ Red ■ Black ■ Blue White

1. Varies with user weight, terrain type, battery amp hour (AH), battery charge, battery condition and tyre condition. This specification can be subject to a variance of 10%.
2. Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of 3%.
3. Tested in accordance with ANSI/RESNA, WC Vol 2, Section 4 and ISO 7176-4 standards. Results derived from the theoretical calculations based on battery specification and drive system performance. Test conducted at maximum weight capacity.
4. Does not account for anti-tip clearance.
5. Width of base (not including seat).
6. Includes standard seat weight.
7. Does not account for anti-tip clearance.
8. Solid durable plastic.
9. Deep cycle (AGM or gel cell type recommended).

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of 3%. Special order seats or non-standard dimensions will not fold-flat and depth-adjustable.