



<21 km
Range per Charge



204 kg
Weight Capacity



55.9 cm
Turning Radius



Jazzy®
EVO 614HD

Heavy-Duty Performance!

The Jazzy® EVO 614HD is designed to offer an unparalleled level of freedom and versatility, empowering users to effortlessly navigate through tight indoor spaces while also conquering outdoor terrains with ease. The Jazzy® EVO 614HD empowers you to pick up and go. That's because it comes engineered with the capabilities that push boundaries of everyday life for the heavy-duty power wheelchair class, with a 204 kg weight capacity.





Heavy Duty Seat Post



Optional Swing Away Controller Mount



Large Footplate with Removable Cover

RANGE PER CHARGE (MAX) ^{6,4}	Up to 21 km
MAXIMUM WEIGHT CAPACITY	204 kg
TURNING RADIUS ¹	55.9 cm
MAXIMUM SPEED ²	7.8 km/h
LENGTH & WIDTH ¹	L: 101.6 cm W: 64.5 cm
WEIGHT OF HEAVIEST PIECE ³	65 kg
GROUND CLEARANCE ¹	7.6 cm at mid-frame
SEAT SIZE WIDTH & DEPTH (STANDARD)	50.8 cm x 50.8 cm (20" x 20" - 22") or 55.8 cm x 55.8 cm (22" x 20" - 22")
STANDARD 18" X 18" SEAT WEIGHT	23.5 kg (high-back)
HIGH-BACK SEAT-TO-FLOOR RANGE	55.8 cm - 60.9 cm
SPECIALTY SEATING OPTIONS	N/A

STANDARD ELECTRONICS	60A, Dynamic Linx
SPECIALTY CONTROLS	N/A
DRIVE TRAIN	Two motor, in-line HD, mid-wheel drive
BRAKING SYSTEM	Regenerative & electro-mechanical
SUSPENSION TYPE	Active-Trac® (ATX)
FRONT WHEELS	6" castors
DRIVE WHEELS	13.75" solid
REAR WHEELS	6" castors
BATTERY REQUIREMENTS	(2) 55AH, 12V
BATTERY CHARGER	Off-board, 5A
BATTERY WEIGHT ⁵	16.8 kg each
COLOUR OPTIONS	<div style="display: flex; align-items: center; gap: 10px;"> <div style="width: 15px; height: 15px; background-color: black; border: 1px solid black;"></div> Black, <div style="width: 15px; height: 15px; background-color: #00AEEF; border: 1px solid black;"></div> Robins Egg Blue, <div style="width: 15px; height: 15px; background-color: #C00000; border: 1px solid black;"></div> Cardinal Red </div>

1. Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.
2. Varies with user weight, terrain, battery condition and tire pressure.
3. Does not include seat or battery weights.
4. Tested in accordance with ANSI/RESNA, WC Vol. 2, Section 4 and ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Testing conducted at maximum weight capacities. This specification can be subject to a variance of +10%, - 5%.
5. Battery weight may vary based on manufacturer.
6. Varies with user weight, terrain, battery type, battery condition and tire pressure. *International configuration shown.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of 3%. Special order seats or non-standard dimensions will not fold-flat and depth-adjustable.