



**<18 km**  
Range per Charge



**181 kg**  
Weight Capacity



**132.5 cm**  
Turning Radius



**PATHRIDER® 10 2.0**

## Built for Style, Comfort and Performance

The new and improved Pathrider® 10 2.0 brings a sleek, refreshed design with a high-back seat and front modesty panel for added style and comfort. Built for stability and power, it delivers everything you'd expect from Pride Mobility. Its four 10-inch pneumatic tyres easily handle varied terrain, ensuring a smooth ride wherever you go. With an adjustable swivel seat, LED headlight, easy-drive tiller, and front basket, the Pathrider® 10 2.0 is the perfect companion for any adventure.





USB Charging Port & Storage Tray



LED Lighting Package



Modesty Panel

RANGE PER CHARGE (MAX) <sup>1,3</sup>	Up to 18 km
MAXIMUM WEIGHT CAPACITY	181 kg
TURNING RADIUS <sup>2</sup>	132.5 cm
MAXIMUM SPEED <sup>1</sup>	8.4 km/h
LENGTH & WIDTH <sup>2,5</sup>	L: 120 cm W: 56.5 cm
SEAT TYPE WIDTH & DEPTH (STANDARD)	High-back seat <sup>8</sup> W: 43.2 cm (17") D: 43.2 cm (17")
SEAT-TO-DECK HEIGHT <sup>2</sup>	43.8 cm - 47.6 cm
SEAT-TO-GROUND HEIGHT <sup>2</sup>	60 cm - 63.8 cm
GROUND CLEARANCE <sup>2,7</sup>	6 cm at motor

OVERALL WEIGHT (NOT INC BATTERY) <sup>2,6</sup>	63 kg
WEIGHT OF HEAVIEST PIECE <sup>2</sup>	29.5 kg (front section)
SUSPENSION	N/A
3W/4W OPTION	4W
FRONT TYRES	10" pneumatic
REAR TYRES	10" pneumatic
BATTERY REQUIREMENTS	(2) 31AH
BATTERY CHARGER	Off-board, 3.5A (32AH)
BATTERY WEIGHT <sup>4</sup>	10.3 kg
COLOUR OPTIONS	■ Red

1. Varies with user weight, terrain type, battery amp hour (AH), battery charge, battery condition and tyre condition. This specification can be subject to a variance of 10%.
2. Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of 3%.
3. Tested in accordance with ANSI/RESNA, WC Vol 2, Section 4 and ISO 7176-4 standards. Results derived from the theoretical calculations based on battery specification and drive system performance. Test conducted at maximum weight capacity. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of 3%.
4. Battery weight may vary based on manufacturer.
5. Width of base (not including seat).
6. Includes standard seat weight.
7. Does not account for anti-tip clearance.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of 3%. Special order seats or non-standard dimensions will not fold-flat and depth-adjustable.