



<25 km Range per Charge



159 kg

Weight Capacity



125 cmTurning Radius



Embrace the Wild!

The Pride Outback 4x4 Powerchair is your ultimate adventure companion, ready to tackle almost any terrain with confidence. Designed for superior performance, this rugged powerchair features powerful 4-pole sealed motors and all-wheel drive for unbeatable traction and stability. Its robust construction ensures a smooth and comfortable ride on uneven surfaces, featuring advanced (SRS) suspension and off-road pneumatic 14" knobby tyres. Whether you're navigating gravel paths, bush tracks, or country properties, the Outback 4x4 Powerchair makes it a breeze.











LED Lighting Package



Centre Mount or Swing Away Leg Rests

RANGE PER CHARGE (MAX) ^{6,4}	Up to 25.75 km at 136 kg Up to 27.4 km at 99.8 kg
MAXIMUM WEIGHT CAPACITY	159 kg
TURNING RADIUS ^{1,7}	125 cm
MAXIMUM SPEED ²	10 km/h
LENGTH ⁷ & WIDTH ¹	L: 106.7 cm W: 73.2 cm
WEIGHT OF HEAVIEST PIECE ³	76.5 kg
GROUND CLEARANCE ¹	10.7 cm at motors
SEAT SIZE WIDTH & DEPTH (STANDARD)	46.5 cm x 46.5 cm (18" x 18")
HIGH-BACK SEAT-TO-FLOOR RANGE	60 cm

1.	Due to manufacturing tolerances and continual product improvement, this
	specification can be subject to a variance of + or - 3 %. Dimensions listed
	are for power base only Overall measurements will vary based on seating and accessory selections.
	and accessory selections.

- 2. Varies with user weight, terrain, battery condition and tire pressure.
- Does not include seat or battery weights.
- Tested in accordance with ANSI/RESNA, WC Vol. 2, Section 4 and ISO 7176-4 standards. Results derived from theoretical calculations based on battery specifications and drive system performance. Testing conducted at maximum weight capacities. This specification can be subject to a variance of +10%, 5%.
- 5. Battery weight may vary based on manufacturer.
- Varies with user weight, terrain, battery type, battery condition and tyre pressure. *International configuration shown.
- 7. Without front riggings

STANDARD ELECTRONICS	120A, Q-Logic 3e
DRIVE TRAIN	4 Motor, independent drive
BRAKING SYSTEM	Regenerative & electro-mechanical
SUSPENSION TYPE	Automotive-grade smooth ride suspension
DRIVE WHEELS	14" knobby pneumatic (38.1 cm)
BATTERY REQUIREMENTS	(2) 95AH, 12V
BATTERY CHARGER	Off-board, 8A
BATTERY WEIGHT⁵	26.6 kg each
COLOUR OPTIONS	Sage Tan

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of 3%. Special order seats or non-standard dimensions will not fold-flat and depth-adjustable.



