



<16 km
Range per Charge



136 kg
Weight Capacity



136.6 cm
Turning Radius



GO GO®
CARBON

Exceptional Portability

The Go Go® Carbon is the perfect blend of strength and portability. Built with a lightweight carbon fibre frame, this 4-wheel scooter offers exceptional durability while weighing just 16 kg without its (lithium) battery. It folds down in seconds for convenient transport and storage, and can be towed like luggage through airports and other environments. Powered by an NFC reader and digital display screen the Go Go® Carbon is designed for a lifestyle on the go.





Conveniently Folds to be Towed like Luggage and Easily Stored



Vibrant Colour Display on Console with NFC Ignition for Extra Security



Light Weight Carbon Fiber Frame

RANGE PER CHARGE (MAX) ^{1,3}	Up to 16 km
MAXIMUM WEIGHT CAPACITY	136 kg
TURNING RADIUS ²	134.6 cm
MAXIMUM SPEED ¹	6 km/h
LENGTH & WIDTH ²	L: 102.3 cm W: 44 cm
SEAT TYPE WIDTH & DEPTH (STANDARD)	Compact, foldable W: 35.6 cm (14") D: 27.5 cm (11")
SEAT-TO-DECK HEIGHT ²	43 cm
SEAT-TO-GROUND HEIGHT ²	57 cm
GROUND CLEARANCE ²	3.8 cm at motor

OVERALL WEIGHT (NOT INC. BATTERY) ²	17.7 kg
WEIGHT OF HEAVIEST PIECE ²	15.9 kg
SUSPENSION	N/A
3W/4W OPTION	4W
FRONT TYRES ²	7" solid
REAR TYRES ²	8" solid
BATTERY REQUIREMENTS	24V 12 Ah lithium (288Wh)
BATTERY CHARGER	Off-board, 2A lithium-ion
BATTERY WEIGHT ²	1.84 kg
COLOUR OPTIONS	■ Black ■ Yellow

1. Varies with user weight, terrain type, battery amp hour (AH), battery charge, battery condition and tyre condition. This specification can be subject to a variance of 10%.
2. Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of 3%.
3. Results derived from theoretical calculations based on lithium battery specifications and drive system performance. Test conducted at maximum weight capacity.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of 3%. Special order seats or non-standard dimensions will not fold-flat and depth-adjustable.



Transport of this item on commercial aircraft is subject to federal regulations and airline policy. Please consult your airline prior to your flight to ensure your Go Go[®] Carbon and its lithium-ion battery are permitted.