

# Airline Travel and Battery Information



## Power Chairs

### Jazzy Carbon Ultra-lite

10.4AH (252Wh) Lithium-ion (standard)  
20AH (504Wh) Lithium-ion (optional)

**Airline Compliant**  
**Not Compliant**



### Jazzy Passport

18.2AH, 24V Lithium-ion (436.8Wh) (standard)  
10.4AH battery (optional)

**Not Compliant**  
**Airline Compliant**



## Scoters

### Go Go Carbon

(1) 24V 12AH lithium (288Wh)

**Airline Compliant**



### Go Go Endurance Li

16AH lithium (2 x 8AH) (385Wh)

**Airline Compliant**



### Quest

(1) 11.5AH Lithium-ion (241.2Wh)

**Airline Compliant**



Transport of this item on commercial aircraft is subject to federal regulations and airline policy. Please consult your airline prior to your flight to ensure your Pride lithium-ion battery are permitted.