



<16.7 km
Range per Charge



136 kg
Weight Capacity



94 cm
Turning Radius



GO GO
ELITE TRAVELLER® 2

Even More Elite!

The Go Go Elite Traveller 2 has the tightest turning radius of all 4-wheel scooters, measuring just 94 cm. With Pride's patented iTurn Technology, you get the best of both worlds: the stability of a 4-wheel scooter and the maneuverability of a 3-wheel model. This unique feature allows you to move smoothly through your home, the store, and other locations with ease and assurance.



GOGO ELITE TRAVELLER[®] 2



Pride's iTurn Technology



Feather Touch Dissassembly



Non Scuffing Flat Free Tyres

RANGE PER CHARGE (MAX) ^{1,3}	Up to 16.7 km
MAXIMUM WEIGHT CAPACITY	136 kg
TURNING RADIUS ²	94 cm
MAXIMUM SPEED ¹	6.8 km/h
LENGTH & WIDTH ^{2,5}	L: 103.6 cm W: 49.4 cm
SEAT TYPE WIDTH & DEPTH (STANDARD)	Compact, foldable ⁷ W: 43 cm (17") D: 43 cm (17")
SEAT-TO-DECK HEIGHT ²	41.6 cm - 45.3 cm
SEAT-TO-GROUND HEIGHT ²	46.3 cm - 51.4 cm
GROUND CLEARANCE ^{2,4}	3.4 cm at transaxle

OVERALL WEIGHT (NOT INC BATTERY) ^{2,6}	37.6 kg
WEIGHT OF HEAVIEST PIECE ²	15.8 kg (front section)
SUSPENSION	N/A
3W/4W OPTION	4W
FRONT TYRES ²	2" x 7" solid
REAR TYRES ²	2.5" x 8" solid
BATTERY REQUIREMENTS ⁸	(2) 18AH
BATTERY CHARGER	Off-board, 2A
BATTERY WEIGHT ²	6 kg
COLOUR OPTIONS	■ Red ■ Black

1. Varies with user weight, terrain type, battery amp hour (AH), battery charge, battery condition and tyre condition. This specification can be subject to a variance of 10%.
2. Due to manufacturing tolerances and continual product improvement, this specification can be subject to variance of 3%.
3. Tested in accordance with ANSI/RESNA, WC Vol 2, Section 4 and ISO 7176-4 standards. Results derived from the theoretical calculations based on battery specification and drive system performance. Test conducted at maximum weight capacity.
4. Does not account for anti-tip clearance.
5. Width of base (not including seat).
6. Includes standard seat weight.
7. Durable solid plastic.

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tyre pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of 3%. Special order seats or non-standard dimensions will not fold-flat and depth-adjustable.

PRIDE[®]
MOBILITY
LIVE YOUR BEST[®]